



Data Sheet

HUMAN BRAIN NATRIURETIC PEPTIDE [BNP]

(aa 1-10)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1034.1 / AE1034.2
Immunogen:	Synthetic human BNP (aa 1-10) poly Lysine conjugated (SPKMOVQSGSGC(G))
Synonyms:	Brain natriuretic peptide 32, BNP-32
Swiss-Prot No:	P16860
Gene Information:	Gene Name: NPPB Gene ID: 4879
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human BNP (aa 1-10)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:500) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2022

For research use only

Publishing research using AE1034? Please let us know so that we can cite your publication as a reference.

