

## Data Sheet

# HUMAN NT-PRO-BNP (aa 1-76)

## SYNTHETIC

<b>Catalog no.:</b>	A 1001AG.1 / A 1001AG.2
<b>Name:</b>	Human N-terminal Pro-Brain Natriuretic Peptide (aa 1-76)
<b>Synonyms:</b>	Natriuretic peptides B, Gamma-brain natriuretic peptide
<b>Swiss-Prot No:</b>	P16860
<b>Gene Information:</b>	Gene Name: NPPB GeneID: 4879
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95 % (HPLC)
<b>Contents:</b>	10 µg / 100 µg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	4 May 2022

**For research use only**

**Publishing research using A 1001AG? Please let us know so that we can cite your publication as a reference.**

