

Data Sheet

HUMAN PRO-BRAIN NATRIURETIC PEPTIDE [BNP] (aa 32-57)

SYNTHETIC

Catalog no.:	AE1027AG.1 / AE1027AG.2
Name:	Human proBNP (aa 32-57)
Synonyms:	Natriuretic peptides B, Gamma-brain natriuretic peptide
Sequence:	QVEQTSLEPLQESPRPTGVWKSREVA
MW:	2.95 kDA
Swiss-Prot No:	P16860
Gene Information:	Gene Name: NPPB GeneID: 4879
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg /1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	27 February 2024

For research use only

Publishing research using AE1027AG? Please let us know so that we can cite your publication as a reference.

