



Data Sheet

HUMAN C-PEPTIDE

ANTIBODY, POLYCLONAL

Catalog no.:	A 5510.1 / A 5510.2
Immunogen:	Synthetic human C-Peptide (EAEDLQVGQVELGGGPGAGSLQPLALEGSLQ) coupled to bovine thyroglobulin
Synonyms:	Insulin (aa 57-87)
Swiss-Prot No:	P01308
Gene Information:	Gene Name: INS GeneID: 3630
Host:	Rabbit
Matrix:	Serum
Specificity	Human and murine C-Peptide
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:2000), ELISA (1:8000), Western Blot (1:100), immunohistochemistry (paraffin sections, 1:2000; cryosections, 1:100), immunocytochemistry (1:100-200) Retrieve antigen by heating paraffin sections in citrate buffer pH 6.0.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	6 April 2022

For research use only

Publishing research using A5510? Please let us know so that we can cite your publication as a reference.

