



Data Sheet

C-PEPTIDE

ANTIBODY, MONOCLONAL

Catalog no.:	A 5513.1 / A 5513.2
Immunogen:	Human C-peptide
Synonyms:	Insulin (aa 57-87)
Swiss-Prot No:	P01308
Gene Information:	Gene Name: INS GenelD: 3630
Host:	Mouse Balb/c
Isotype:	IgG ₁
Matrix:	50 mM Na-citrate pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃
Specificity:	Human C-peptide Kd = 7x 10 ⁻⁹ mol/l Cross reactivity with human proinsulin was 0.93%.
Contents:	10 µg / 100 µg (liquid)
Known applications:	ELISA (capture) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (do not freeze)
Last updated on:	21 February 2024

For research use only

Publishing research using A 5513? Please let us know so that we can cite your publication as a reference.

