



## Data Sheet

# HUMAN PRO C-TYPE NATRIURETIC PEPTIDE [proCNP] (aa 51-79) ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE 1002.1 / AE 1002.2
<b>Immunogen:</b>	Synthetic human proCNP (aa 51-79) (DLRVDTKSRAAWARLLQEHPNARKYKGAN)
<b>Swiss-Prot No:</b>	P23582
<b>Gene Information:</b>	Gene Name: NPPC, CNP2 GenelD: 4880
<b>Host:</b>	Sheep
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Synthetic human proCNP (aa 51-79)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:1000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	14 April 2022

**For research use only**

**Publishing research using AE1002? Please let us know so that we can cite your publication as a reference.**

