



Data Sheet

HUMAN CALCITONIN

ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | A 4111.1 / A 4111.2 |
| Immunogen: | Synthetic human Calcitonin, BSA-conjugated (CGNLSTCMLGTYTQDFNKFHTFPQTAIGVGAP) |
| Swiss-Prot No: | P01258 |
| Gene Information: | Gene Name: CALCA, CALC1 GeneID: 796 |
| Host: | Sheep |
| Matrix: | Serum |
| Specificity: | Human Calcitonin There was no cross reactivity obtained with human Katalcalcin or human Calcitonin Gene related Peptide. |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | RIA (1:80 000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 18 March 2022 |

For research use only

Publishing research using A 4111? Please let us know so that we can cite your publication as a reference.

