



Data Sheet

HUMAN CALCITONIN GENE-RELATED PEPTIDE 2

ANTIBODY, POLYCLONAL

Catalog no.:	A 4612.1 / A 4612.2
Immunogen:	Synthetic human Calcitonin gene-related peptide 2, bTG-conjugated (ACNTATCVTHRLAGLLSRSGGMVKSNTFVPTNVGSKAF)
Synonyms:	Beta-type CGRP, Beta-CGRP, Calcitonin gene-related peptide II, CGRP-II
Swiss-Prot No:	P10092
Gene Information:	Gene Name: CALCB, CALC2 GeneID: 797
Host:	Goat
Matrix:	Serum
Specificity:	Human Calcitonin Gene-related Peptide There was no cross reactivity obtained with human and salmon Katalcalcin or Calcitonin.
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:3000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	6 April 2022

For research use only

Publishing research using A 4612? Please let us know so that we can cite your publication as a reference.

