

## Data Sheet

# HUMAN CALCITONIN GENE-RELATED PEPTIDE 1 SYNTHETIC

<b>Catalog No.:</b>	A 4613AG.1 / A 4613AG.2
<b>Name:</b>	Synthetic human Calcitonin gene-related peptide 1
<b>Synonyms:</b>	Alpha-type CGRP, Calcitonin gene-related peptide I, CGRP-I
<b>Sequence:</b>	ACDTATCVTHRLAGLLSRSGGVKNNFVPTNVGSKAF
<b>MW:</b>	3.79 kDa
<b>Swiss-Prot No:</b>	P06881
<b>Gene Information:</b>	Gene Name: CALCA, CALC1 Gene ID: 796
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration.
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	05 May 2022

**For research use only**

**Publishing research using A 4613AG? Please let us know so that we can cite your publication as a reference.**

