

Data Sheet

HUMAN MRP8/14 (CALPROTECTIN, L1 PROTEIN) RECOMBINANT

Catalog no.:	A 6940AG.1 / A 6940AG.2
Name:	Recombinant human Calprotectin, tetrameric heterocomplex of MRP8 and MRP14
Synonyms:	Calprotectin, L1 protein, S100A8/A9 complex, calgranulin A/calgranulin B complex
Swiss-Prot No:	P05109 and P06702
Gene Information:	Gene Name: S100A8 (MRP8) and S100A9 (MRP14) Gene ID: 6279 and 6280
Source:	<i>E. coli</i>
Storage buffer:	Proprietary; contains < 0.1 % NaN ₃
Purity:	> 90% (capillary gel electrophoresis)
Known applications:	ELISA This antigen has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Contents:	10 µg / 100 µg (liquid)
Store at:	2-8 °C (do not freeze)
Last updated on:	24 May 2022

For research use only

Publishing research using A 6940AG? Please let us know so that we can cite your publication as a reference.

