



Data Sheet

COW MILK CASEIN

ANTIBODY, POLYCLONAL

Catalog no.:	A 1220.3 / A 1220.4
Immunogen:	Native casein purified from cow milk
Host:	Goat
Matrix:	IgG, Protein G purified from serum; PBS pH 7,4; 0,095 % NaN ₃
Specificity:	Cow milk casein
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:5000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	17 March 2022

For research use only

Publishing research using A 1220? Please let us know so that we can cite your publication as a reference.

