



## Data Sheet

# HUMAN CHEMOKINE CC-1

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE 1119.1 / AE 1119.2
<b>Immunogen:</b>	Native human Chemokine CC-1, poly Lysin conjugated
<b>Synonyms:</b>	None
<b>Swiss-Prot No:</b>	Q76BR7
<b>Gene Information:</b>	Gene Name: CCL1 GeneID: 6346
<b>Host:</b>	Chicken
<b>Matrix:</b>	IgY
<b>Specificity:</b>	Human Chemokine CC-1
<b>Contents:</b>	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	12 January 2018

### For research use only

**Publishing research using AE1119? Please let us know so that we can cite your publication as a reference.**

