



## Data Sheet

# HUMAN CHEMOKINE CC-2 (aa 78-92) SYNTHETIC

|                          |                                                                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog no.:</b>      | AE1115AG.1 / AE1115AG.2                                                                                                                                   |
| <b>Name:</b>             | Human Chemokine CC-2 (aa 78-92)                                                                                                                           |
| <b>Synonyms:</b>         | HCC-2, C-C motif chemokine 15, Leukotactin-1 (LKN-1), MIP-1 delta, Macrophage inflammatory protein 5 (MIP-5), Mrp-2b, NCC-3, Small-inducible cytokine A15 |
| <b>Sequence:</b>         | PGVQDCMKKLKPYSI                                                                                                                                           |
| <b>MW:</b>               | 1.71 kDa                                                                                                                                                  |
| <b>Swiss-Prot No:</b>    | Q16663                                                                                                                                                    |
| <b>Gene Information:</b> | Gene Name: CCL15<br>Gene ID: 6359                                                                                                                         |
| <b>Source:</b>           | Synthetic                                                                                                                                                 |
| <b>Purity:</b>           | > 95% (HPLC)                                                                                                                                              |
| <b>Contents:</b>         | 100 µg / 1 mg (lyophilized)<br>Resuspend 0.1% acetic acid for required concentration.                                                                     |
| <b>Store at:</b>         | - 20 °C<br>Repeated thawing and freezing must be avoided                                                                                                  |
| <b>Last updated on:</b>  | 2 April 2025                                                                                                                                              |

**For research use only**

**Publishing research using AE1115AG? Please let us know so that we can cite your publication as a reference.**

