



Data Sheet

ANTI HUMAN CHOLECYSTOKININ-RECEPTOR [CCK-R] TYPE A (aa 197-210)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog no.: | AE1130.1 / AE1130.2 |
| Immunogen: | Synthetic human CCK-R Type A (aa 197-210) poly Lysin conjugated (RFLPPNDVMQQSWH) |
| Synonyms: | CCK-A receptor, CCK-AR, Cholecystokinin-1 receptor, CCK1-R |
| Swiss-Prot No: | P32238 |
| Gene Information: | Gene Name: CCKAR, CCKRA GenelD: 886 |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Human CCK-R Type A (aa 197-210) |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | Western Blot (1:2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 22 April 2022 |

For research use only

Publishing research using AE1130? Please let us know so that we can cite your publication as a reference.

