

## Data Sheet

# HUMAN CHOLECYSTOKININ-RECEPTOR TYPE A [CCK-R] (aa 180-195)

## SYNTHETIC

<b>Catalog no.:</b>	AE1129AG.1 / AE1129AG.2
<b>Name:</b>	Human CCK-R Type A (aa 180-195)
<b>Synonyms:</b>	CCK-A receptor (CCK-AR), Cholecystokinin-1 receptor (CCK1-R)
<b>Sequence:</b>	SNLVPFTKNNNQTANM
<b>MW:</b>	1.79 kDa
<b>Swiss-Prot No:</b>	P32238
<b>Gene Information:</b>	Gene Name: CCKAR, CCKRA Gene ID: 886
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	19 May 2022

**For research use only**

**Publishing research using AE1129AG? Please let us know so that we can cite your publication as a reference.**

