



## Data Sheet

# HORSE CHORIONIC GONADOTROPIN (CG)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AS1004.1 / AS1004.2
<b>Immunogen:</b>	Purified horse Chorionic Gonadotropin
<b>Synonyms:</b>	Pregnant Mare Serum Gonadotropin (PMSG)
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Horse Chorionic Gonadotropin
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	RIA (1:500.000), Western Blot (1:500-10.000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	19 September 2024

### For research use only

**Publishing research using AS1004? Please let us know so that we can cite your publication as a reference.**

