



## Data Sheet

# HUMAN CHROMOGRANIN B [CgB] (aa 480-493)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AST1020.1 / AST1020.2
<b>Immunogen:</b>	Synthetic human CgB (aa 480-493) KLH-conjugated ( <u>C</u> RFQDKQYSSHTAE) Underlined amino acids were substituted or added for technical reasons.
<b>Synonyms:</b>	Secretogranin-1, Chromogranin-B, CgB, Secretogranin I, Sgl
<b>Swiss-Prot No:</b>	P05060
<b>Gene Information:</b>	Gene Name: CHGB GeneID: 1114
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human CgB (aa 480-493)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	RIA, Western Blot, immunohistochemistry (paraffin sections)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	4 May 2022

### For research use only

**Publishing research using AST1020? Please let us know so that we can cite your publication as a reference.**

