



Data Sheet

HUMAN DERMCIDIN PEPTIDE (aa 66-80)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1065.1 / AE1065.2
Immunogen:	Synthetic human Dermcidin Peptide (aa 66-80), polyLysine-conjugated (LEKGLDGAKKAVGGL)
Synonyms:	Preproteolysin
Swiss-Prot No:	P81605
Gene Information:	Gene Name: DCD GeneID: 117159
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Dermcidin Peptide (aa 66-80)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:30 000), RIA (1:5000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE1065? Please let us know so that we can cite your publication as a reference.

