

Data Sheet

HUMAN DERMICIDIN PEPTIDE (aa 96-110)

SYNTHETIC

| | |
|--------------------------|--|
| Catalog no.: | AE1061AG.1 / AE1061AG.2 |
| Name: | Synthetic human Dermcidin Peptide (aa 96-110), polyLysine-conjugated Underlined amino acids differ from the original sequence; they were added or substituted for technical reasons |
| Synonyms: | Preproteolysin |
| Sequence: | KGAVHDVKDVLDSVL, KGAVHDVKDVLDSVL <u>gK7A</u> |
| MW: | 13600 Da |
| Swiss-Prot No: | P81605 |
| Gene Information: | Gene Name: DCD GeneID: 117159 |
| Source: | Synthetic |
| Contents: | 100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration |
| Store at: | - 20 °C Repeated thawing and freezing must be avoided |
| Last updated on: | 18 May 2022 |

For research use only

Publishing research using AE1061AG? Please let us know so that we can cite your publication as a reference.

