



Data Sheet

HUMAN DIAMINE OXIDASE (DAO)

ANTIBODY, MONOCLONAL

Catalog no.:	A 7003.1 / A 7003.2
Immunogen:	Recombinant human Diamine oxidase (DAO)
Synonyms:	Amiloride-sensitive amine oxidase (copper-containing), Amiloride-binding protein (ABP), Histaminase Kidney amine oxidase (KAO), EC 1.4.3.22
Swiss-Prot No:	P19801
Gene Information:	Gene Name: ABP1, AOC1, DAO1 GenelD: 26
Host:	Mouse Balb/c
Clone no.:	FP96
Isotype:	IgG ₁
Matrix:	Cell culture supernatant, Protein G purified, PBS pH 7.4
Specificity:	Diamine oxidase
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	7 April 2022

For research use only

Publishing research using A 7003? Please let us know so that we can cite your publication as a reference.

