



Data Sheet

HUMAN DIAMINE OXIDASE (DAO)

ANTIBODY, POLYCLONAL

Catalog no.:	A 7001.1 / A 7001.2
Immunogen:	Recombinant human Diamine oxidase (DAO)
Synonyms:	Amiloride-sensitive amine oxidase (copper-containing), Amiloride-binding protein (ABP), Histaminase Kidney amine oxidase (KAO), EC 1.4.3.22
Swiss-Prot No:	P19801
Gene Information:	Gene Name: ABP1, AOC1, DAO1 GeneID: 26
Host:	Rabbit
Matrix:	Serum
Specificity:	Diamine oxidase
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:4000), Western Blot (1:2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	7 April 2022

For research use only

Publishing research using A 7001? Please let us know so that we can cite your publication as a reference.

