


Data Sheet

HUMAN DIAMINE OXIDASE (DAO)

RECOMBINANT

Catalog no.:	A 7003AG.1 / A 7003AG.2
Name:	Recombinant human Diamine oxidase (DAO)
Synonyms:	Amiloride-sensitive amine oxidase (copper-containing), Amiloride-binding protein (ABP), Histaminase Kidney amine oxidase (KAO), EC 1.4.3.22
MW:	83.3 kDa
Swiss-Prot No:	P19801
Gene Information:	Gene Name: ABP1, AOC1, DAO1 GenelD: 26
Source:	Produced in eukaryotic cell line
Purification:	Heparin sepharose affinity chromatography
Storage buffer:	PBS pH 7.2
Contents:	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest. for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	06 May 2022



For research use only

Publishing research using A 7003AG? Please let us know so that we can cite your publication as a reference.



Immundiagnostik AG

Stubenwald-Allee 8a · 64625 Bensheim · Germany

Phone: +49 6251 70190-0 · Fax: +49 6251 70190-363 · dept.immuochemicals@immundiagnostik.com · www.immundiagnostik.com