



## Data Sheet

# HUMAN ENDOTHELIN-1

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AS1015.1 / AS1015.2
<b>Immunogen:</b>	Human Endothelin-1
<b>Swiss-Prot No:</b>	P05305
<b>Gene Information:</b>	Gene name: EDN1 GeneID: 1906
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Endothelin-1, pig and domestic ruminant Endothelin-1
<b>Contents:</b>	20µl / 100 µl (lyophilized) Resuspend in aqua bidest. for required concentration
<b>Known applications:</b>	RIA, EIA, Western Blot, immunohistochemistry (paraffin sections)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	28 April 2022

### For research use only

**Publishing research using AS1015? Please let us know so that we can cite your publication as a reference.**

