



Data Sheet

HUMAN FIBRONECTIN

NATIVE

Catalog no.:	A 4290AG.1 / A 4290AG.2
Name:	Fibronectin, FN
Synonyms:	Cold-insoluble globulin
MW:	221 kDa
Swiss-Prot No:	P02751
Gene Information:	Gene FN1 GeneID: 2335
Source:	Human Plasma
Purification:	Affinity chromatography
Storage buffer:	0.1 M Tris pH 7.4
Contents:	250 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	05 May 2022

For research use only

Publishing research using A 4290AG? Please let us know so that we can cite your publication as a reference.

