



## Data Sheet

# HUMAN FIBRONECTIN, CELLULAR

## NATIVE

<b>Catalog no.:</b>	A 4291AG.1 / A 4291AG.2
<b>Name:</b>	Fibronectin, FN
<b>Synonyms:</b>	Cold-insoluble globulin
<b>MW:</b>	560 kDa
<b>Swiss-Prot No:</b>	P02751
<b>Gene Information:</b>	Gene FN1 GeneID: 2335
<b>Source:</b>	Human foreskin fibroblasts
<b>Purification:</b>	Affinity chromatography
<b>Storage buffer:</b>	PBS pH 7.4
<b>Contents:</b>	250 µg / 1 mg (lyophilized) Resuspend in 250 µl / 1 ml aqua bidest at 37 °C
<b>Known applications:</b>	Cell culture substrate
<b>Store at:</b>	- 80 °C ??? (31.05.2021 MKA) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	5 May 2022

### For research use only

**Publishing research using A 4291AG? Please let us know so that we can cite your publication as a reference.**

