



## Data Sheet

# HUMAN GLUCAGON-LIKE PEPTIDE-2 [GLP-2] SYNTHETIC

<b>Catalog No.:</b>	A 6102AG.1 / A 6102AG.2
<b>Name:</b>	Human GLP-2
<b>Sequence:</b>	HADGFSDEMNTILDNLAARDFINWLIQTKITD
<b>Swiss-Prot No:</b>	P01275
<b>Gene Information:</b>	Gene Name: GCG GeneID: 2641
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 90% (HPLC)
<b>Contents:</b>	100 µg / 1mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	06 May 2022

**For research use only**

**Publishing research using A 6102AG? Please let us know so that we can cite your publication as a reference.**

