



Data Sheet

RAT GUANYLIN (aa 78-92)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1019.1 / AE1019.2
Immunogen:	Synthetic rat Guanylin (aa 78-92) active peptide (PNTCEICAYAACTGC)
Synonyms:	Guanylate cyclase activator 2A
Swiss-Prot No:	P28902
Gene Information:	Gene Name: Guca2a, Guca2 GeneID: 25656
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic rat Guanylin (aa 78-92) Recognizes the intracellular form of Guanylin (aa 1-115) and the native circulating peptide in rat blood
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:3000); RIA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	20 April 2022

For research use only

Publishing research using AE1019? Please let us know so that we can cite your publication as a reference.

