



Data Sheet

HUMAN PRO-GUANYLIN (aa 13-25)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1114.1 / AE1114.2
Immunogen:	Synthetic human pro-Guanylin (aa 13-25) (ESVKKLKDLEPEQ)
Synonyms:	Guanylate cyclase activator 2A (aa 13-25), Guanylate cyclase-activating protein I (aa 13-25), GCAP-I (aa 13-25)
Swiss-Prot No:	Q02747
Gene Information:	Gene Name: GUCA2A, GUCA2 GenelD: 2980
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human pro-Guanylin (aa 13-25)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA, RIA, Western Blot This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE1114? Please let us know so that we can cite your publication as a reference.

