

Data Sheet

HUMAN pro-GUANYLIN [aa 13-25]

SYNTHETIC

Catalog no.:	AE1114AG.1 / AE1114AG.2
Name:	Human pro-Guanylin (aa 13-25)
Synonyms:	Guanylate cyclase activator 2A (aa 13-25), Guanylate cyclase-activating protein I (aa 13-25), GCAP-I (aa 13-25)
Sequence:	ESVKKLKLQEPQ
Swiss-Prot No:	Q02747
Gene Information:	Gene Name: GUCA2A, GUCA 2 GeneID: 2980
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	19 May 2022

For research use only

Publishing research using AE1114AG? Please let us know so that we can cite your publication as a reference.

