

Data Sheet

HUMAN pro-GUANYLIN [aa 78-94]

SYNTHETIC

Catalog no.:	AE1009AG.1 / AE1009AG.2
Name:	Human pro-Guanylin (aa 78-94)
Sequence:	EDPGTCEICAYAACTGC
Swiss-Prot No:	Q02747
Gene Information:	Gene Name: GUCA2A, GUCA 2 GeneID: 2980
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	12 May 2022

For research use only

Publishing research using AE1009AG? Please let us know so that we can cite your publication as a reference.

