

Data Sheet

RAT GUANYLIN (aa 78-92)

SYNTHETIC

Catalog no.:	AE1019AG.1 / AE1019AG.2
Name:	Rat Guanylin (aa 78-92) active peptide
Synonyms:	Guanylate cyclase activator 2A
Sequence:	PNTCEICAYAACTGC
Swiss-Prot No:	P28902
Gene Information:	Gene Name: Guca2a, Guca2 GenelD: 25656
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	12 May 2022

For research use only

Publishing research using AE1019AG? Please let us know so that we can cite your publication as a reference.

