



Data Sheet

HEMOGLOBIN

ANTIBODY, POLYCLONAL

Catalog no.:	A 7817.1 / A 7817.2
Immunogen:	Human Hemoglobin
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Hemoglobin
Contents:	20 µl / 100 µl Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

Last updated on: 7 April 2022

For research use only

Publishing research using A 7817? Please let us know so that we can cite your publication as a reference.

