



## Data Sheet

# HUMAN HEMOGLOBIN ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 7818.1 / A 7818.2
<b>Immunogen:</b>	Human Hemoglobin
<b>Host:</b>	Sheep
<b>Matrix:</b>	Immunoglobulins, antigen-affinity purified from serum; PBS
<b>Specificity:</b>	Human Hemoglobin
<b>Contents:</b>	10 µg / 100 µg Resuspend in 10 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided.
<b>Last updated on:</b>	7 April 2022

**For research use only**

**Publishing research using A 7818? Please let us know so that we can cite your publication as a reference.**

