



Data Sheet

INFLIXIMAB

ANTIBODY, RECOMBINANT MONOCLONAL, COUPLED TO HSA

Catalog no.:	A 1230.3 / A 1230.4
Immunogen:	Infliximab
Synonyms:	Remicade®, Inflectra®, Remsima®
CAS-No.:	170277-31-3
Host:	Human cell line
Clone no.:	AF312-A05-HSA
Form:	Monovalent single chain variable fragment (scFv) coupled to HSA
Matrix:	PBS, 0.09% NaN ₃
Specificity:	Infliximab variable domain
Contents:	10 µg / 100 µg (liquid)
Known applications:	ELISA (capture)
	<small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C
Last updated on:	31 October 2024

For research use only

Publishing research using A 1230? Please let us know so that we can cite your publication as a reference.

