



Data Sheet

HUMAN INSULIN-LIKE GROWTH FACTOR-2 [IGF-2]

ANTIBODY, POLYCLONAL

Catalog no.:	AS1009.1 / AS1009.2
Immunogen:	Human Insulin-like Growth Factor-2 (Somatomedin A)
Synonyms:	Somatomedin-A, T3M-11-derived growth factor
Swiss-Prot No:	P01344
Gene Information:	Gene Name: IGF2 GenelD: 3481
Host:	Rabbit
Matrix:	Serum
Specificity:	Human and domestic ruminants IGF-2
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	Western Blot (1:1000-2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	28 April 2022

For research use only

Publishing research using AS1009? Please let us know so that we can cite your publication as a reference.

