



## Data Sheet

# HUMAN INSULIN-LIKE GROWTH FACTOR-BINDING PROTEIN 2 [IGFBP-2] (aa 168-185)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE1060.1 / AE1060.2
<b>Immunogen:</b>	Synthetic human IGFBP-2 (aa 168-185), polyLysine-conjugated (GKGGKHHHLGLEEPKLRP)
<b>Swiss-Prot No:</b>	P18065
<b>Gene Information:</b>	Gene Name: IGFBP2, BP2, IBP2 GenelD: 3485
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human IGFBP-2 (aa 168-185)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	22 April 2022

### For research use only

**Publishing research using AE1060? Please let us know so that we can cite your publication as a reference.**

