



Data Sheet

BOVINE INTERFERON α 1 (IFN α 1)

ANTIBODY, POLYCLONAL

Catalog no.:	AS1020.1/AS1020.2
Immunogen:	Recombinant bovine Interferon α 1
Host:	Rabbit
Matrix:	Serum
Specificity:	Bovine Interferon α 1
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	Western Blot This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	28 April 2022

For research use only

Publishing research using AS1020? Please let us know so that we can cite your publication as a reference.

