

Data Sheet

HUMAN INTERLEUKIN-1 β [IL-1 β]

ANTIBODY, POLYCLONAL

Catalog no.:	AS1018.1 / AS1018.2
Immunogen:	Recombinant human Interleukin-1 β
Synonyms:	IL-1 beta, Catabolin
Swiss-Prot No:	P01584
Gene Information:	Gene Name: IL1B, IL1F2 GeneID: 3553
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Interleukin-1 β
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	Western Blot <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	19 September 2024

For research use only

Publishing research using AS1018? Please let us know so that we can cite your publication as a reference.

