



Data Sheet

BOVINE INTERLEUKIN-2 (IL-2)

ANTIBODY, POLYCLONAL

Catalog no.:	AS1023.1 / AS1023.2
Immunogen:	Recombinant bovine Interleukin-2
Synonyms:	T-cell growth factor (TCGF)
Swiss-Prot No:	P05016
Gene Information:	Gene Name: IL2, IL-2 GeneID: 280822
Host:	Rabbit
Matrix:	Serum
Specificity:	Bovine Interleukin-2
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	Western Blot This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	28 April 2022

For research use only

Publishing research using AS1023? Please let us know so that we can cite your publication as a reference.

