



Data Sheet

RAT GUANYLATE CYCLASE 2G (aa 1045-1057)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1125.1 / AE1125.2
Immunogen:	Synthetic rat Guanylate Cyclase 2G (aa 1045-1057) polyLysine-conjugated (EFTEEEAKVPEIL)
Synonyms:	Guanylyl cyclase receptor G (GC-G), Kinase-like domain-containing soluble guanylyl cyclase (ksGC)
Swiss-Prot No:	P55205
Gene Information:	Gene Name: Gucy2g, Ksgc GeneID:
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic rat Guanylate Cyclase 2G (aa 1045-1057)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	Western Blot (1:1000-2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE1125? Please let us know so that we can cite your publication as a reference.

