



Data Sheet

HUMAN ALPHA KLOTHO

ANTIBODY, RECOMBINANT MONOCLONAL

Catalog no.:	A 1101.1 / A 1101.2
Immunogen:	Human secreted α Klotho (isoform 2, aa 1-516)
Synonyms:	EC: 3.2.1.31
Swiss-Prot No:	Q9UEF7
Gene Information:	Gene Name: KL Gene ID: 9365
Host:	Human cell line
Clone no.:	E10
Form:	Monovalent single-chain variable fragment fused to mouse IgG2a constant domain (scFv-IgG2a)
Matrix:	PBS pH 7.4, 0.09% NaN ₃
Specificity:	Human Klotho
Contents:	10 μ g / 100 μ g (liquid)
Known applications:	ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	23 August 2022

For research use only

Publishing research using A 1101? Please let us know so that we can cite your publication as a reference.

