



## Data Sheet

# HUMAN LUTEINIZING HORMONE [LH]

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AS1000.1/AS1000.2
<b>Immunogen:</b>	Bovine Luteinizing Hormone
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Luteinizing Hormone
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	RIA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved)
<b>Last updated on:</b>	27 April 2022

### For research use only

**Publishing research using AS1000? Please let us know so that we can cite your publication as a reference.**

