



## Data Sheet

# HUMAN LYSOZYME

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 6905.1 / A 6905.2
<b>Immunogen:</b>	Human lysozyme
<b>Synonyms:</b>	Lysozyme C, renal amyloidosis
<b>Swiss-Prot No:</b>	P61626
<b>Gene Information:</b>	Gene Name: LYZ GenelD: 4069
<b>Host:</b>	Sheep
<b>Matrix:</b>	IgG, affinity purified using human lysozyme; PBS
<b>Specificity</b>	Human lysozyme
<b>Contents:</b>	10µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	7 April 2022

### For research use only

**Publishing research using A 6905? Please let us know so that we can cite your publication as a reference.**

