



## Data Sheet

# MISTLETOE LECTIN

## ANTIBODY, MONOCLONAL

<b>Catalog no.:</b>	A 1087.1 / A 1087.2
<b>Immunogen:</b>	Recombinant mistletoe lectin
<b>Host:</b>	Mouse Balb/c
<b>Clone no.:</b>	1C5-H8
<b>Isotype:</b>	IgM
<b>Matrix:</b>	Cell culture supernatant, ThioSorb M purified; PBS pH 7.4
<b>Specificity:</b>	Recombinant mistletoe lectin
<b>Contents:</b>	10 µg / 100 µg Resuspend in 10 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	16 March 2022

### For research use only

**Publishing research using A 1087? Please let us know so that we can cite your publication as a reference.**

