



Data Sheet

MISTLETOE LECTIN

ANTIBODY, MONOCLONAL

Catalog no.:	A 1086.1 / A 1086.2
Immunogen:	Recombinant mistletoe lectin
Host:	Mouse Balb/c
Clone no.:	2G8-F7
Isotype:	IgM
Matrix:	Cell culture supernatant, ThioSorb M purified; PBS pH 7.4
Specificity:	Recombinant mistletoe lectin
Contents:	10 µg / 100 µg Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1086? Please let us know so that we can cite your publication as a reference.

