



## Data Sheet

# HUMAN MYELOPEROXIDASE

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 6840.1 / A 6840.2
<b>Immunogen:</b>	Human Myeloperoxidase
<b>Synonyms:</b>	MPO
<b>Swiss-Prot No:</b>	P05164
<b>Gene Information:</b>	Gene Name: MPO GenelD: 4353
<b>Host:</b>	Rabbit
<b>Matrix:</b>	IgG, Protein G purified from serum; TRIS pH 7.4
<b>Specificity:</b>	Human Myeloperoxidase
<b>Contents:</b>	10 µg / 100 µg (lyophilized) Resuspend in 10 µ / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA, immunohistochemistry (paraffin sections) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.</small>
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	06 April 2022

### For research use only

**Publishing research using A 6840? Please let us know so that we can cite your publication as a reference.**

