

Data Sheet

HUMAN PRO-MYOSTATIN (aa 79-92)

ANTIBODY, POLYCLONAL

Catalog no.:	AK3004.1 / AK3004.2
Immunogen:	Synthetic human Pro-Myostatin (aa 79-92) KLH conjugated (SDGSLEDDDYHATT)
Synonyms:	Growth/differentiation factor 8 (GDF-8)
Swiss-Prot No:	O14793
Gene Information:	Gene Name: MSTN, GDF8 Gene ID: 2660
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Pro-Myostatin (aa 79-92)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:500) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	27 April 2022

For research use only

Publishing research using AK3004? Please let us know so that we can cite your publication as a reference.

