



## Data Sheet

# HUMAN MYOSTATIN

## ANTIBODY, POLYCLONAL, BIOTINYLATED

<b>Catalog no.:</b>	AK3009.5 / AK3009.6
<b>Immunogen:</b>	Recombinant human Myostatin
<b>Synonyms:</b>	Growth / differentiation factor 8 (GDF-8)
<b>Swiss-Prot No:</b>	O14793
<b>Gene Information:</b>	Gene Name: MSTN, GDF8 GeneID: 2660
<b>Host:</b>	Goat
<b>Matrix:</b>	IgG, Protein G purified from serum; PBS, 0.09% NaN <sub>3</sub>
<b>Specificity:</b>	Human Myostatin
<b>Contents:</b>	10 µg / 100 µg (liquid)
<b>Known applications:</b>	ELISA
	This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C
<b>Last updated on:</b>	27 April 2022

**For research use only**

**Publishing research using AK3009? Please let us know so that we can cite your publication as a reference.**

