



Data Sheet

HUMAN MYOSTATIN

ANTIBODY, POLYCLONAL

Catalog no.:	AK3009.3 / AK3009.4
Immunogen:	Recombinant human Myostatin
Synonyms:	Growth / differentiation factor 8 (GDF-8)
Swiss-Prot No:	O14793
Gene Information:	Gene Name: MSTN, GDF8 GeneID: 2660
Host:	Goat
Matrix:	IgG, Protein G purified from serum; 200 mM TRIS, pH 8.0
Specificity:	Human Myostatin
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	27 April 2022

For research use only

Publishing research using AK3009? Please let us know so that we can cite your publication as a reference.

