

Data Sheet

HUMAN MYOSTATIN (aa 83-98)

SYNTHETIC

Catalog no.:	AK 3005AG.1 / AK 3005AG.2
Name:	Human Myostatin (aa 83-98)
Sequence:	NMLYFNGKEQIIYGKI
Synonyms:	Growth/differentiation factor 8, GDF-8
Swiss-Prot No:	O14793
Gene Information:	Gene name: MSTN, GDF8 GenelD: 2660
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	20 May 2022

For research use only

Publishing research using AK3005AG? Please let us know so that we can cite your publication as a reference.

